


#### **KEY FEATURES**

- Complex wavefront modulation
   63 electrodes enabling replication of up to the 7th radial order of Zernike polynomials (>35 modes)
   with high fidelity
- Straightforward system integration
   Compact housing compatible with standard 30 mm cage systems by rods, lens tubes, and post assemblies
- Linear & hysteresis-free response
   Electrostatic actuation suited for open-loop wavefront control
- Remarkable optical quality
   Active best flat with an induced
   RMS wavefront error of less than λ/40
- Polarization-independent
   Wavefront modulation independent of the light polarization for maximized efficiency







## **TECHNOLOGY**

The Delta 7 is based on the Deformable Phase Plate (DPP) technology, exclusively developed by Phaseform GmbH. DPP is composed of a fluidic chamber, enclosed by a thin membrane, which is deformed by electrostatic force. The force is generated by a 2D array of transparent electrodes embedded within the optical aperture of the DPP. The sophisticated optofluidic design of the DPP enables gravity-neutral performance for orientation-independent, high-quality wavefront modulation.

## **GENERAL SPECIFICATIONS**

Modulator type

Clear aperture diameter

Number of actuators

Number of actuators across aperture diameter

Connectivity

Operating system

**Driving software** 

Included in the Delta 7 package

Optofluidic DPP (Deformable Phase Plate), electrostatically actuated

10 mm

63

7

**USB 2.0** 

Windows, Linux, and macOS

SDK and GUI available in Python. Wrapper to execute Python functions in Matlab.

Driving electronics, control software, cables, manual

## **OPTICAL SPECIFICATIONS**

Wavefront RMS error of best flat

Maximum peak-to-valley of the generated wavefronts

Maximum spatial frequency of the correction

Optical transmission (VIS-NIR version)

Laser Induced Damage Threshold (LIDT)

Nominal operation laser power

< 15 nm (orientation independent)

 $> 8 \mu m$ 

7th radial order of Zernike modes

425 nm - 2200 nm 75% at  $\lambda$ =800 nm (no AR coatings applied)

10 W/cm<sup>2</sup> for 10s @ 1070nm CW

Factory calibrated for < 100 mW CW (over full optical aperture)



## **MECHANICAL SPECIFICATIONS**

Thickness (within clear aperture)

Response time (best flat to maximum deformation)

Hysteresis

Linearity

Mounting compatibility

Connector cable length

0.87 mm

 $< 55 \, \mathrm{ms}$ 


<1%


> 92 %

30 mm cage system rods, SM1 tubing, and ø=1/2" post

1.5 m, optional: 0.3 m

#### HOUSING MECHANICAL DRAWINGS





(New version coming soon...)

## **ELECTRICAL SPECIFICATIONS**

**Actuator voltage** 

Maximum power consumption

Power supply

up to 300 V DC

< 9 W

120/230 VAC, 2.5 phono plug (included)

## THERMAL SPECIFICATIONS

Storage temperature

Operating temperature

10 °C to 35 °C


20 °C to 25 °C




# PEAK TO VALLEY OF THE GENERATED ZERNIKE MODES (OPTICAL PATH DIFFERENCES IN µm)

Z ( 1,-1)	8.0	Z (4,-4)	2.0	Z (5,3)	1.1	Z (7,-7)	0.7
Z ( 1,1)	8.0	Z (4,-2)	1.8	Z (5,5)	1.3	Z (7,-5)	0.6
Z (2,-2)	5.1	Z (4,0)	1.5	Z (6,-6)	1.0	Z (7,-3)	0.5
Z (2 ,0)	5.0	Z (4,2)	1.8	Z (6,-4)	0.9	Z (7,-1)	0.5
Z (2,2)	5.1	Z (4,4)	2.0	Z (6,-2)	0.9	Z (7,1)	0.5
Z (3,-3)	3.6	Z (5,-5)	1.3	Z (6,0)	0.7	Z (7,3)	0.5
Z (3,-1)	3.0	Z (5,-3)	1.1	Z (6,2)	0.9	Z (7,5)	0.6
Z (3,1)	3.0	Z (5,-1)	1.1	Z (6,4)	0.9	Z (7,7)	0.7
Z (3,3)	3.6	Z (5,1)	1.1	Z (6,6)	1.0		
Z (3,1)	3.0	Z (5,-1)	1.1	Z (6,4)	0.9	• • •	

#### **RESPONSE TIME**







## DISCLAIMER

10%

20

40

60

Time [ms]

0

0

(Simulation model in OpticStudio ZEMAX is available)

All specifications are preliminary and subject to change without further notice. No representation or warranty, either expressed or implied, is made as to the reliability, completeness, or accuracy of this specification sheet.

100

80

## **FURTHER INFORMATION**

Phaseform GmbH Georges-Köhler-Allee 302 D-79110 Freiburg i.B., Germany

Email: info@phaseform.com Web: www.phaseform.com Phone: +49 761 6007 9018

For more information visit our website and send us an email to request a live demonstration.