

# DELTA 7<sup>10</sup>

## TRANSMISSIVE WAVEFRONT MODULATOR

### DPP TECHNOLOGY

The Delta 7 is based on the Deformable Phase Plate (DPP) technology, exclusively developed by Phaseform GmbH. DPP is composed of a fluidic chamber, enclosed by a thin membrane, which is deformed by electrostatic force. The force is generated by a 2D array of transparent electrodes embedded within the optical aperture of the DPP. The sophisticated optofluidic design of the 10 mm aperture DPP enables gravity-neutral performance for orientation-independent, high-quality wavefront modulation.

### KEY FEATURES

#### **Complex wavefront modulation**

63 electrodes enabling replication of up to the 7th radial order of Zernike polynomials (>35 modes) with high fidelity

#### **Straightforward system integration**

Compact housing compatible with standard 30 mm cage systems by rods, lens tubes, and post assemblies

#### **Linear & hysteresis-free response**

Electrostatic actuation suited for open-loop wavefront control

#### **Remarkable optical quality**

Active best flat with an induced RMS wavefront error of less than  $\lambda/40$  independent of the DPP orientation

#### **Polarization-independent**

Wavefront modulation independent of the light polarization for maximized efficiency



# SPECIFICATIONS

## GENERAL

Modulator type

Optofluidic DPP (Deformable Phase Plate),  
electrostatically actuated

Clear aperture diameter

10 mm

Number of actuators

63

Number of actuators across aperture diameter

7

Connectivity

USB 2.0

Included in the Delta 7 package

Driving electronics, control software, cables, manual

## OPTICAL

Wavefront RMS error of best flat

< 15 nm (orientation independent)

Maximum peak-to-valley of the generated wavefronts

> 7  $\mu\text{m}$

Maximum spatial frequency of the correction

7th radial order of Zernike modes

Wavefront RMS drift

< 5% after 60 min

Laser Induced Damage Threshold (LIDT)


10 W/cm<sup>2</sup> for 10s @ 1070 nm CW

Optical transmission

Without AR coating




With AR coated membrane



## GENERATED ZERNIKE MODES IN OPEN LOOP

Max amplitude RMS, peak-to-valley, and purity of generated Zernike modes

Z (n,m)	RMS [ $\mu\text{m}$ ]	PV [ $\mu\text{m}$ ]	Purity [%]	Z (n,m)	RMS [ $\mu\text{m}$ ]	PV [ $\mu\text{m}$ ]	Purity [%]
Z (1,-1)	2	7	98	Z (5,3)	0.15	0.7	82
Z (1,1)	2	7	97	Z (5,5)	0.15	0.9	90
Z (2,-2)	1	3.8	96	Z (6,-6)	0.1	0.8	83
Z (2,0)	1	3.8	96	Z (6,-4)	0.07	0.5	70
Z (2,2)	1	3.8	96	Z (6,-2)	0.07	0.5	77
Z (3,-3)	0.5	2.5	94	Z (6,0)	0.1	0.6	77
Z (3,-1)	0.5	1.7	94	Z (6,2)	0.07	0.5	73
Z (3,1)	0.5	1.7	93	Z (6,4)	0.07	0.5	70
Z (3,3)	0.5	2.5	95	Z (6,6)	0.1	0.8	77
Z (4,-4)	0.25	1.3	93	Z (7,-7)	0.07	0.6	74
Z (4,-2)	0.25	1.1	88	Z (7,-5)	0.07	0.5	70
Z (4,0)	0.25	1.2	91	Z (7,-3)	0.07	0.5	70
Z (4,2)	0.25	1.1	86	Z (7,-1)	0.07	0.5	70
Z (4,4)	0.25	1.3	94	Z (7,1)	0.07	0.5	70
Z (5,-5)	0.15	0.9	90	Z (7,3)	0.07	0.5	70
Z (5,-3)	0.15	0.7	84	Z (7,5)	0.07	0.5	70
Z (5,-1)	0.15	0.8	80	Z (7,7)	0.07	0.6	78
Z (5,1)	0.15	0.8	80				



\*Purity is defined as the fraction of the target Zernike mode relative to the root-square-sum of all modes.

\*\*The Zernike-tree shows experimentally measured replicated Zernike modes up to the 7th order, in open-loop, using an interferometric setup.

# SPECIFICATIONS, CONT.

## MECHANICAL

Thickness (within clear aperture)  
 Response time (best flat to maximum deformation)  
 Hysteresis  
 Linearity  
 Mounting orientation

0.87 mm  
 < 40 ms  
 < 1%  
 > 92%  
 Arbitrary (horizontal and vertical standard, other orientations upon request)  
 30 mm cage system rods, SM1 tubing, and  $\varnothing=1/2"$  post  
 1.5 m

Mounting capability  
 Connector cable length

## ELECTRICAL

Actuator voltage  
 Maximum power consumption  
 Power supply

up to 295 V DC  
 < 9 W  
 120/230 VAC

## THERMAL

Storage temperature  
 Operating temperature

10 °C to 35 °C  
 20 °C to 25 °C

# SOFTWARE

## Graphical User Interface (GUI)




- Python-based interface using serial communication
- Direct application of phase profiles using Zernike polynomial descriptions

## Software Development Kit (SDK)

- Python scripting control interface
- MATLAB wrappers to execute Python functions
- Easy integration into custom control software
- Reference implementation of the serial communication protocol included

## Phinden (optional add-on for Delta 7)


- Turnkey software for real-time, image-based adaptive optics in microscopy
- Automated aberration correction improving resolution, SNR, and contrast
- Runs seamlessly alongside native microscope control software



## System Requirements

Operating system: Windows, Linux, macOS  
 Python: Minimum version 3.8, recommended - 3.11  
 MATLAB: Minimum version R2020b


# OPTICS HOUSING MECHANICAL DRAWINGS



## APPLICATION IDEA

Using a Delta 7<sup>10</sup> in the excitation path of a light-sheet microscope to expand the effective field of view via modulating the Gaussian beam profile, and a Delta 7<sup>20</sup> in the detection path for wavefront coding to extend depth of field.

*Courtesy of Prof. Pablo Loza-Alvarez and Dr. Gustavo Castro, SLN Lab, ICFO, Spain*



## DISCLAIMER

All specifications are preliminary and subject to change without notice. No representation or warranty, either expressed or implied, is made as to the reliability, completeness, or accuracy of this specification sheet.

## CONTACT US

Get in touch with us today for a live demo!

**Phaseform GmbH**  
Georges-Köhler-Allee 302  
79110 Freiburg i.B.  
Germany

[www.phaseform.com](http://www.phaseform.com)  
[info@phaseform.com](mailto:info@phaseform.com)  
[+49 761 216 0800 0](tel:+4976121608000)



Phaseform is supported by

European  
Innovation  
Council

